Sorbitol/Mannitol 3% B. Braun

bg	Инструкции за употреба	4
da	Brugsanvisning	6
de	Gebrauchsanweisung	8
en	Information for use	10
es	Instrucciones de uso	12
fi	Käyttöohje	14
fr	Mode d'emploi	16
hr	Upute za uporabu	18
it	Istruzioni per l'uso	20
lt	Vartojimo instrukcija	22
mk	Упатство за користење	24
nl	Gebruiksaanwijzing	26
no	Brukerveiledning	28
pt	Instruções de uso	30
sr	Uputstvo za upotrebu	32
sv	Bruksanvisning	34

Black

Dimension: 148 x 210 mm
38 pages

LLD-Spec.: L48

Production code: ---

Medical Device___Irrig...sung___0193_GIF
0193/NU19399/1224
GIF (A5-Booklet, 16-spr.)
Production site: Crissier


Font size: 7,0 pt.
6,0 pt. (Table)

V-0100


B | BRAUN

Device description


Ecobag® CE standard:


Ecobag® CE click:

Port system Ecobag® CE standard:

Port system Ecobag® CE click:

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This packaging unit of sale contains multiple containers of irrigation solution. Such containers shall not be distributed individually unless an instruction for use is provided with each container.

Ecobag® CE standard/ Ecobag® CE click: Sorbitol/Mannitol 3% B. Braun is stored in an irrigation bag and is applied with the help of an irrigation set.

- Ecobag® CE standard: Irrigation bag with two identical access ports. Ecobag® CE standard port can be used with irrigation sets.
- Ecobag® CE standard port cannot be used with Ecoclick irrigation sets.
- Ecobag® CE click: Irrigation bag with a spike-free connector to the irrigation set. Ecobag® CE click port must be used with spike-free B. Braun Ecoclick® irrigation sets.

Sorbitol/Mannitol 3% B. Braun is an irrigation medium intended to be used for intra- and postoperative irrigation purposes. Sorbitol and Mannitol are constituents of the irrigation solution to adjust the osmolarity of the solution in order to reduce haemolysis, but without being electrically conductive which would interfere with high frequency (HF) endoscopic surgery. Do not re-use.

Sterility

Sterile and endotoxin-free irrigation solution. Sterilization is performed by moist heat.

Do not use if package is damaged.

Visually inspect the sterile barrier system [for breaches and integrity prior to use](#).

Intended purpose

Sorbitol/Mannitol 3% B. Braun is a ready for use, sterile, single use medical device, intended to be used for irrigation applications (irrigation solution).

Indication

Sorbitol/Mannitol 3% B. Braun is used as an irrigation solution in patients undergoing:

- endoscopic examinations of body cavities
- urological interventions by endoscopy involving HF (high frequency) current, e.g. TUR (transurethral resections), disintegration of bladder stones, interventions of the urether, interventions of the renal pelvis etc.
- gynaecological interventions by endoscopy involving HF current
- arthroscopy of the knee, shoulder or other joints involving HF current
- postoperative irrigation

Patient population

Sorbitol/Mannitol 3% B. Braun can be used for all patients under consideration of the respective indications, contraindications and warnings for use.

Safety and efficacy in pediatric patients have not been established.

Safety and efficacy in pregnant women and geriatric patients have not been established.

Intended User

To be used by country-specific authorized persons e.g.:

Healthcare professionals (professional nurses, licensed and/or certified physicians, healthcare workers, midwives, paramedics, pharmacists).

Contraindications

Do not use for intravascular infusion, injection and/or intrathecal irrigation.

Sorbitol/Mannitol 3 % B. Braun must not be applied in patients with:

- oligo-anuria
- fructose-sorbitol intolerance
- hypersensitivity to any of the substances or materials employed

Residual Risks/Side effects

Side effects:

Side effects are possible.

Sorbitol/Mannitol 3% B. Braun contains sorbitol. If there is undiagnosed hereditary fructose-sorbitol intolerance the administration of solutions containing sorbitol can lead to nausea, hypoglycaemia, elevation of lactate levels, acidosis and acute liver damage, with possible lethal consequences. The most commonly observed side effect is the TURP (transurethral resection of the prostate) syndrome (the absorption of fluid with hyponatraemia), caused by an excessive inflow of fluid into the circulatory system.

This may result in:

- Renal insufficiency (especially where there is existing renal damage)
- Haemolysis
- Clotting disorders caused by inhibition of the clotting factor activity or by lowering the coagulation factor concentration through haemodilution
- Hypothermia

TURP syndrome is characterized by the following symptoms:

- Nausea, vomiting
- Neurological problems: headache, sleepiness, agitation, confusion, which can develop into a potentially fatal coma or convulsions
- Visual disturbances
- Hyponatraemia, hypocalcaemia, sometimes pronounced, reduced osmotic pressure, reduced levels of haematocrit and serum protein
- Bradycardia
- Increased blood volume indicated by dyspnoea, variations in blood pressure or even acute pulmonary or cerebral oedema

All these symptoms are linked to an excessive inflow of irrigation liquid into the systemic circulation and indicate hyperhydration which can be fatal. They can appear during intervention (such as long-lasting interventions, particularly when a large venous sinus is opened) or in the decremental phase.

Often seen after transurethral resection of the prostate (TURP), the described symptoms and signs have been referred to as the TURP syndrome, although they have also been described after other urological or gynaecological surgical procedures.

If this diagnosis is suspected, these symptoms require the intervention to be stopped after careful haemostasis. The naatraemia and the haematocrit should be rapidly measured and appropriate treatment for hyperhydration with haemodilution should begin – notably by restricting the hydric flow and administering a suitable diuretic.

Precautions

It is necessary to exclude the possibility of a hereditary fructose intolerance or a fructose-1,6-diphosphatase deficiency before commencement of irrigation. If details of the case history are not available or are not unequivocal a fructose intolerance test must be carried out first. If this is not done there is the possibility of severe complications (including a lethal outcome under certain circumstances).

Only use the solution if the tamper evidence seal/container is undamaged and the solution is clear and colourless. Otherwise discard container.

Plastic containers can occasionally be damaged during transport from the manufacturer to the hospital, but also within the clinic. This can lead to contamination with microbial or fungiform growth of the solution. A careful visual inspection is therefore always required before using the solution. Particular attention must be paid even to minimal damage of the sterility barrier (i.e. to the closure) or missing integrity of the tamper evidence seal due to possible risk of contamination.

Sorbitol/Mannitol 3% B. Braun is not intended to be combined or used with medicinal products. Healthcare professional facing use of Sorbitol/Mannitol 3% B. Braun irrigation solution in conjunction with medicinal products must check the safety information of the additive as provided by the respective manufacturer.

The container can be warmed up at 40°C for maximum 2 weeks, and must not be warmed above 40°C. Temperature of the product must be checked before usage based on the clinical need.

The equipment used for warming must allow to regulate the required temperature.

Warning

Sorbitol/Mannitol 3 % B. Braun should be applied with care in:

- cardiopulmonary diseases
- disturbances of renal function.

In order to prevent TURP syndrome, the following recommendations should be followed:

- Time of surgery should be reduced to the minimum especially for those procedures at the increased risk of excessive inflow (normally no excess of 60 minutes)
- Pressure in the operating cavity should be kept at a minimum
- Ensure adequate haemostasis
- Sufficient monitoring of fluid balance and electrolytes

Throughout the procedure, frequent and regular checks must be carried out for excessive inflow of irrigation liquid into the systemic circulation. This is particularly suggested by the appearance of any of the symptoms characteristic for the TURP syndrome during local or regional anaesthetics (for detail please refer to section side effects). The symptoms can also appear during the decremental phase of an intervention.

Disturbances of the acid-base balance and disturbances of the water electrolyte balance should be corrected before the use of Sorbitol/Mannitol 3 % B. Braun.

Strictly sterile precautions are to be taken when handling the irrigation solution.

Sorbitol/Mannitol 3% B. Braun is intended for single use only and does not have to be completely emptied.

Opened containers must not be saved, quantities left over must be discarded.

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability.

Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Store out of reach children.


The solution must not be used after the expiry date given on the label.

Operating Instructions

Ecobag® CE standard:

1. Remove the protective cap of the port of Ecobag® CE Irrigation Container.
2. Remove the protective cap of the Irrigation Set and close the roller clamp.
3. Pre-insertion: Introduce the spike in the center of the gum, perpendicularly to the surface, until you feel a resistance. Do not insert completely the spike to avoid wedging the gum in the port.
4. Slight Rotation and relaxation :
 - Slightly rotate the spike until you feel that the gum has broken.
 - Relax the strength given on the spike. The gum must come back at its original position.
5. Final insertion : Introduce completely the spike into the port.
6. Visual inspection : Before using the solution, always inspect the content of the bag for visible particles.


Ecobag CE standard

Ecobag CE click

Ecobag® CE click :

1. Remove the protective cap of the port of Ecobag® CE click Irrigation Container.
2. Remove the protective cap of the Irrigation Set and close the roller clamp.
3. Connect the port of Ecobag® CE click Irrigation container with the Ecoclick® Irrigation Set.

Duration of use

The volumes used and the duration of use depend on the extent of the intervention and on the postoperative course.

In operative settings, the continuous treatment with the solutions should not exceed 60 minutes to prevent TURP syndrome.

Disposal

Disposal according to local guidelines and / or clinical protocols.

Storage and handling conditions

Keep away from sunlight.

Keep dry (unless it is to be warmed up in a water bath immediately before use).

Temperature (T) range between 5°C and 25°C.

The following disinfectants are compatible with the outer surface of Ecobag :

- 2-Propanol 70%
- 80% Ethanol
- Chlorhexidine 2% (e.g. Softacloth®)
- Ethanol 74% + Propan-2-ol 10% (e.g. Softasept® N)
- 1-Propanol 70%
- Sterile 70% denatured ethanol blended with deionised water
- Mixture of isopropyl alcohol with hydrogen peroxide
- H₂O₂ 30%

Clinical benefits & Performances

The summary of safety and clinical performance of the product is available in the European database on medical devices (EUDAMED) at ec.europa.eu/tools/eudamed using the following Basic UDI-DI:

Sorbitol/Mannitol 3 B. Braun Irrigation Ecobag (40392390000024522L)

The clinical benefit of the irrigation for the patient depends on many variables including the mechanism of fluid delivery, site and type of intervention, wound type, patient condition, and is generally directly related to the result of the therapy. Therefore, therapy effect is an indirect patient benefit of irrigation solution when used as intended.

Notice to the user

All serious incidents related to this product must be reported to B. Braun and the competent authority in the country where the product is being used.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

Date of issue and if appl. revision date

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B. Braun Medical AG
B. Braun Melsungen AG