ENGLISH	FRANÇAIS	DEUTSCH	ITALIANO	ESPAÑOL	PORTUGUÊS DO BRASIL	NEDERLANDS	DANSK
 smith&nephew	 smith&nephew	 smith&nephew	 smith&nephew	 smith&nephew	 smith&nephew	 smith&nephew	 smith&nephew

RAPID RHINO® SINU-KNIT®

Instructions for Use

DESCRIPTION

RAPID RHINO® SINU-KNIT® Dressings are well knitted carbonylmethylcellulose (CMC) hydro fibres. The dressing is an inverted multi layer design that allows the fabric to act as a stent. The hydrocolloid fabric creates a moist gel barrier between the incision and the mucosal tissue due to its hygroscopic capability.

INDICATIONS FOR USE

An intranasal splint is intended to minimize bleeding and edema and to prevent adhesions between the septum and the nasal cavity. It is placed in the nasal cavity after surgery or trauma.

CONTRAINDICATIONS

There are no known contraindications. Medical judgment is necessary when using RAPID RHINO SINU-KNIT Dressings.

⚠️ WARNINGS AND PRECAUTIONS

Do not attempt to use the device for anything other than the indicated use.

Adverse Reactions

- These include the possibility of infection, recurrent or persistent bleeding, and dislodgement of the packing necessitating medical treatment.

- Toxic shock syndrome has been reported to be associated with the use of certain nasal packing materials with fluid absorbing properties. This condition can cause death, hence all patients should be monitored for any signs, which indicate the development of this condition.

- Beautifying devices is a possibility in physically small patients or patients with a shallow nasopharynx due to improper positioning of the device.

- Respiratory distress or decreased oxygenation can occur in patients at risk for impaired respiratory condition.

DIRECTIONS FOR USE

- Remove RAPID RHINO SINU-KNIT Dressing from the sterile packaging and place on a sterile surgical field.
- Grasp the fabric with dry non-sterreded forceps and gently advance the dressing into the desired location in a dry (un-hydrated) state. To ensure the maximum efficacy, be sure to cover the denuded tissue and mucosa.

- Wet the fabric with Sterile Water (not Saline). Ensure that the entire fabric from the deepest portion of the posterior cavity, to the anterior has been totally hydrated allowing the fabric conversion into its gelled state.

- Have the patient irrigate with light saline solution until first office visit (usually 7-10 days post op) Do not irrigate with hypertonic solution until after first office visit.

- Any residual gel that has not dissolved or left the surgical site through postoperative passages may be easily removed through gentle suction.

HOW SUPPLIED

The device is supplied sterile, and is for **SINGLE USE ONLY. Do not clean, resterilize, or reuse the device, as this may damage or compromise the performance resulting in product malfunction, failure, or patient injury.** Cleaning, reesterilization, or reuse of the device may also expose the patient to the risk of transmitting infectious diseases.

The product is sterilized using ethylene oxide.

The product is not made with natural rubber latex and is DEHP free.

DISPOSAL

The device is to be disposed of in accordance with applicable government regulations, and the Medical Waste Management Plan for the facility.

CUSTOMER SERVICE

This product is guaranteed for materials, function, and workmanship for single patient use only. **DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.**

Product Complaints
All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized representative.

	ArthroCare Corporation 7000 West William Cannon Drive Austin, TX 78753 USA Phone: +1 800-343-5717 CUSTOMER SERVICE www.arthrocare.com
	Smith & Nephew Orthopaedics GmbH Altenrainstraße 14 76532 Tullington Germany ec_rep@smith-nephew.com

CE/REP Smith & Nephew Orthopaedics GmbH
Altenrainstraße 14
76532 Tullington
Germany
ec_rep@smith-nephew.com

SMITH
This product is guaranteed for materials, function, and workmanship for single patient use only. **DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.**

Product Complaints
All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized representative.

	Catologue number
	Lot number
	Use-by Date
	Consult Instructions for Use
	Do not reuse
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Keep away from sunlight
	Caution
	Sterilized using ethylene oxide
	Manufacturer
	Authorized Representative in the European Community
	CAUTION: Federal law (USA) restricts this device to sale only by or on the order of a physician.
	CE mark and identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).

*Trademark of Smith & Nephew.
This product may be covered by one or more US Patents. See smith-nephew.com/patents for details.

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© Smith & Nephew 2022.