DERMABOND ADVANCED™

TOPICAL SKIN ADHESIVE

(2-Octyl Cyanoacrylate)

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(2-Octyl Cyanoacrylate)

Description

DERMABOND ADVANCED™ Topical Skin Adhesive is a sterile, liquid topical skin adhesive containing a monomeric (2-octyl cyanoacrylate) formulation and the colorant D&C Violet No. 2. It is provided as a single use applicator in a blister package. The applicator is composed of a crushable glass ampoule contained within a plastic applicator tip. As applied to the skin, the liquid is syrup-like in viscosity and polymerizes within minutes. Studies have shown that following application of DERMABOND ADVANCED™ Topical Skin Adhesive, it acts as a barrier to prevent microbial infiltration into the healing wound. See Directions for use.

DERMABOND ADVANCED™ Topical Skin Adhesive is intended for use only by healthcare professionals who are trained in surgical techniques requiring the use of topical skin adhesives.

Clinical benefits expected from skin approximation are to promote wound healing by primary intention after surgical incisions and simple, thoroughly cleansed, trauma-induced lacerations.

Indications / Intended Use

DERMABOND ADVANCED™ Topical Skin Adhesive is intended for topical application only to hold closed easily approximated skin edges of wounds from surgical incisions, including incisions from minimally invasive surgery, and simple, thoroughly cleansed, trauma-induced lacerations. DERMABOND ADVANCED™ Topical Skin Adhesive may be used in conjunction with, but not in place of, deep dermal stitches.

Contraindications

- Do not use on any wound with evidence of active infection, gangrene, or wounds of decubitus etiology.

- Do not use on mucosal surfaces or across mucocutaneous junctions (e.g., oral cavity, lips), or on skin which may be regularly exposed to body fluids or with dense natural hair (e.g., scalp).
- Do not use on patients with a known hypersensitivity to cyanoacrylate, formaldehyde, or benzalkonium chloride.

Patient Target Group(s)

DERMABOND ADVANCED™ Topical Skin Adhesive is to be used on anyone requiring the use of topical skin adhesives to help approximate skin edges of wound from surgical incisions and simple, thoroughly cleansed, trauma-induced lacerations.

Warnings

- DERMABOND ADVANCED™ Topical Skin Adhesive is a fast setting adhesive capable of adhering to most body tissue and many other materials, such as surgical gloves and stainless steel. Inadvertent contact with any body tissue, and any surfaces or equipment that are not disposable or that cannot be readily cleaned with a solvent such as acetone should be avoided.
- Polymerization of DERMABOND ADVANCED™ Topical Skin Adhesive may be accelerated by water or fluids containing alcohol. DERMABOND ADVANCED™ Topical Skin Adhesive should not be applied to wet wounds.
- DERMABOND ADVANCED™ Topical Skin Adhesive should not be applied to the eye. If contact with the eye occurs, flush the eye copiously with saline or water. If residual adhesive remains, apply topical ophthalmic ointment to help loosen the bond and contact an ophthalmologist.
- When closing facial wounds near the eye with DERMABOND ADVANCED™ Topical Skin Adhesive, position the patient so

that any run-off of adhesive is away from the eye. The eye should be closed and protected with gauze. Prophylactic placement of petroleum jelly around the eye, to act as a mechanical barrier or dam, can be effective in preventing inadvertent flow of adhesive into the eye. DERMABOND ADVANCED™ Topical Skin Adhesive will not adhere to skin pre-coated with petroleum jelly. Therefore, avoid using petroleum jelly on any skin area where DERMABOND ADVANCED™ Topical Skin Adhesive is intended to adhere. Use of DERMABOND ADVANCED™ Topical Skin Adhesive near the eye has inadvertently caused some patients' eyelids to be sealed shut. In some of these cases, general anesthesia and surgical removal has been required to open the eyelid.

- DERMABOND ADVANCED™ Topical Skin Adhesive should not be used below the skin because the polymerized material is not absorbed by tissue and can elicit a foreign body reaction.
- DERMABOND ADVANCED™ Topical Skin Adhesive should not be used in high skin tension areas or across areas of increased skin tension, such as joints e.g., knuckles, elbows, or knees, unless the joint will be immobilized during the skin healing period or unless skin tension has been removed by application of another wound closure device (e.g., sutures or skin staples) prior to application of DERMABOND ADVANCED™ Topical Skin Adhesive.
- DERMABOND ADVANCED™ Topical Skin Adhesive treated wounds should be monitored for signs of infection. Wounds with signs of infection, such as erythema, edema, warmth, pain and pus, should be evaluated and treated according to standard practice for infection.
- DERMABOND ADVANCED™ Topical Skin Adhesive should not be used on wound sites that will be subjected to repeated or prolonged moisture or friction.
- DERMABOND ADVANCED™ Topical Skin Adhesive should only be used after wounds have been adequately and thoroughly cleaned, debrided and internal tissue layers have been closed to relieve skin tension in accordance with standard surgical practice. Local anesthetic may be used when necessary to ensure adequate cleansing and debridement.
- Excessive pressure of the applicator tip against wound edges or surrounding skin can force the wound edges apart and allow adhesive into the wound. Adhesive within the wound could delay wound healing and/or result in adverse cosmetic outcome. Therefore, DERMABOND ADVANCED™ Topical Skin Adhesive should be applied with a very light brushing motion of the applicator tip over easily approximated wound edges.
- DERMABOND ADVANCED™ Topical Skin Adhesive polymerizes through an exothermic reaction in which a small

amount of heat is released. With the proper technique of applying DERMABOND ADVANCED™ Topical Skin Adhesive in one continuous layer onto a dry wound and allowing time for polymerization, heat is released slowly and the sensation of heat or pain experienced by the patient is minimized. However, if DERMABOND ADVANCED™ Topical Skin Adhesive is applied so that large drops of liquid are allowed to remain unspread, the patient may experience a sensation of heat or discomfort.

- DERMABOND ADVANCED™ Topical Skin Adhesive is packaged for single patient use. Discard remaining opened material after each wound closure procedure.
- Do not resterilize/reuse DERMABOND ADVANCED™ Topical Skin Adhesive.
- Do not place DERMABOND ADVANCED™ Topical Skin Adhesive in a procedure pack/tray that is to be sterilized prior to use. Exposure of DERMABOND ADVANCED™ Topical Skin Adhesive to excessive heat (as in autoclaves or ethylene oxide sterilization) or radiation (such as gamma or electron beam), is known to increase its viscosity and may render the product unusable.

Precautions

- Do not apply liquid or ointment medications or other substances to the wound after closure with DERMABOND ADVANCED™ Topical Skin Adhesive, as these substances can weaken the polymerized film and allow for dehiscence (skin edge separation). DERMABOND™ Adhesive permeability by topical medications has not been studied. Prior to application, cleanse the application site thoroughly to remove any remaining blood, fluids, or topical medications/anesthetics.
- DERMABOND™ Adhesive permeability by fluids is not known and has not been studied.
- DERMABOND ADVANCED™ Topical Skin Adhesive, as a liquid, is syrup-like in viscosity. To prevent inadvertent flow of liquid DERMABOND ADVANCED™ Topical Skin Adhesive to unintended areas: (1) the patient should be positioned so that the wound is preferably in a horizontal position, with DERMABOND ADVANCED™ Topical Skin Adhesive applied from above, and (2) DERMABOND ADVANCED™ Topical Skin Adhesive should be applied in one continuous layer rather than in large drops.
- Hold applicator away from yourself and the patient with the tip pointing down and break the ampoule close to the center of the bulb one time only. Do not crush the contents of the applicator repeatedly as further manipulation of the applicator may cause glass shard penetration of the outer tube. Glass shard

penetration can result in inadvertent skin punctures, which may result in the transmission of bloodborne pathogens.

- DERMABOND ADVANCED™ Topical Skin Adhesive should be used immediately after crushing the glass ampoule as the adhesive will polymerize and not flow freely from the applicator, rendering the device unusable.
- If unintended bonding of intact skin occurs, peel the adhesive from the skin, but do not pull the skin edges apart. Petroleum jelly or acetone may help loosen the bond. Other agents such as water, saline, povidone-iodine, chlorhexidine gluconate, or soap, are not expected to immediately loosen the bond.
- Safety and effectiveness of DERMABOND ADVANCED™ Topical Skin Adhesive on wounds of patients with peripheral vascular disease, insulin dependent diabetes mellitus, blood clotting disorders, personal or family history of keloid formation or hypertrophy, burst stellate lacerations, animal or human bites, or puncture or stab wounds, have not been studied.
- Safety and effectiveness on wounds that have been treated with DERMABOND ADVANCED™ Topical Skin Adhesive and then exposed for prolonged periods to direct sunlight or tanning lamps have not been studied.
- Safety and effectiveness of DERMABOND ADVANCED™ Adhesive on wounds on vermilion surfaces has not been studied.

Adverse Reactions / Undesirable Side Effects

Adverse reactions related to either the wound closure procedure or the use of DERMABOND ADVANCED™ Topical Skin Adhesive are possible. The following events have been identified as potentially associated with the wounds closed with DERMABOND ADVANCED™ Topical Skin Adhesive.

- Infection (redness more than 3-5 mm from the wound margin, swelling, purulent discharge, pain, increased skin temperature, fever)
- Acute inflammation (erythema, edema, pain, warmth)
- Dehiscence (Skin Edge Separation)
- Excessive itching
- Skin blistering

Events potentially associated with the wound closure procedure include bleeding, skin edge necrosis, seroma, and hematoma.

Healthcare professionals should convey adverse reactions, undesirable side effects, and risks associated with the product and the procedure to the patient and advise the patient to contact a healthcare professional in case of any deviation from the normal postoperative course. Any serious incident that has occurred in

relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information

DERMABOND ADVANCED™ Topical Skin Adhesive is Magnetic Resonance (MR) Safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogen (H340), mutagen (H350) or reproductive toxicant (H360) based on human evidence and animal studies.

Directions for use

1. The application of DERMABOND ADVANCED™ Topical Skin Adhesive requires thorough wound cleansing. Follow standard surgical practice for wound preparation before application of DERMABOND ADVANCED™ Topical Skin Adhesive (i.e., anesthetize, irrigate, debride, obtain hemostasis and close deep layers).
2. Pat the wound dry with dry, sterile gauze to ensure direct tissue contact for adherence of the DERMABOND ADVANCED™ Topical Skin Adhesive to the skin. Moisture accelerates DERMABOND ADVANCED™ Topical Skin Adhesive's polymerization and may affect wound closure results.
3. To prevent inadvertent flow of liquid DERMABOND ADVANCED™ Topical Skin Adhesive to unintended areas of the body, the patient should be positioned so that the wound is preferably in a horizontal position. DERMABOND ADVANCED™ Topical Skin Adhesive should be applied from above the wound.
4. DERMABOND ADVANCED™ Topical Skin Adhesive should be used immediately after crushing the glass ampoule, since the liquid DERMABOND ADVANCED™ Topical Skin Adhesive will flow freely from the tip for only a few minutes.
5. Refer to the instructions on the blister package for crushing the glass ampoule and expressing the liquid adhesive. Remove the applicator from the blister package. Hold the applicator with the thumb and a finger, away from the patient to prevent any unintentional placement of the liquid DERMABOND ADVANCED™ Topical Skin Adhesive into the wound or on the patient. With the applicator tip pointing downward, apply pressure at the midpoint of the bulb to crush the glass ampoule. Gently squeeze the applicator sufficiently to moisten the internal filter with the liquid adhesive. Stop squeezing and allow the liquid DERMABOND ADVANCED™ Topical Skin Adhesive to draw back into the applicator.

6. Approximate wound edges with gloved fingers or sterile forceps. Slowly apply the liquid DERMABOND ADVANCED™ Topical Skin Adhesive in one continuous layer to the surface of the approximated wound edges using a gentle brushing motion. Maintain manual approximation of the wound edges for approximately 60 seconds after the application. The layer width can be increased or decreased by adjusting the amount of pressure applied to the bulb during application.

NOTE: DERMABOND ADVANCED™ Topical Skin Adhesive polymerizes through an exothermic reaction. If the liquid DERMABOND ADVANCED™ Topical Skin Adhesive is applied so that large drops are allowed to remain without being evenly spread, the patient may experience a sensation of heat or discomfort. The sensation may be higher on sensitive tissues. This can be minimized by applying DERMABOND ADVANCED™ Topical Skin Adhesive in one continuous layer.

NOTE: Excessive pressure of the applicator tip against the wound edges or surrounding skin can result in forcing the wound edges apart and allowing DERMABOND ADVANCED™ Topical Skin Adhesive into the wound. DERMABOND ADVANCED™ Topical Skin Adhesive within the wound could delay wound healing and/or result in adverse cosmetic outcome.

NOTE: Full apposition strength is expected to be achieved within minutes after the adhesive is applied. Full polymerization is expected when the DERMABOND ADVANCED™ Topical Skin Adhesive layer is no longer sticky.

7. Do not apply liquid or ointment medications onto wounds closed with DERMABOND ADVANCED™ Topical Skin Adhesive because these substances can weaken the polymerized film, leading to dehiscence (skin edge separation).
8. Protective dry dressings such as gauze may be applied only after DERMABOND ADVANCED™ Topical Skin Adhesive film is completely solid/polymerized: not tacky to the touch after a few minutes. Allow the DERMABOND ADVANCED™ Topical Skin Adhesive to fully polymerize before applying a bandage. If a dressing, bandage, adhesive backing or tape is applied before complete polymerization, the dressing can adhere to the film. The film can be disrupted from the skin when the dressing is removed, and dehiscence (skin edge separation) can occur.
9. Patients should be instructed to not pick at the polymerized film of DERMABOND ADVANCED™ Topical Skin Adhesive. Picking at the film can disrupt its adhesion to the skin and cause dehiscence (skin edge separation). Picking at the film can be discouraged by an overlying dressing.
10. Apply a dry protective dressing for children or other patients who may not be able to follow instructions for proper wound care.

11. Patients treated with DERMABOND™ Topical Skin Adhesive should be instructed on, how to care for the wound after it's treated with DERMABOND ADVANCED™ Topical Skin Adhesive. This instruction should be reviewed with each patient or guardian to ensure understanding of the proper care for the treatment site.
12. Patients should be instructed that until the polymerized film of DERMABOND ADVANCED™ Topical Skin Adhesive has sloughed naturally (usually in 5-10 days), there should be only transient wetting of the treatment site. Patients may shower and bathe the site gently. The site should not be scrubbed, soaked, or exposed to prolonged wetness until after the film has sloughed naturally and the wound has healed closed. Patients should be instructed not to go swimming during this period.
13. If removal of DERMABOND ADVANCED™ Topical Skin Adhesive is necessary for any reason, carefully apply petroleum jelly or acetone to the DERMABOND ADVANCED™ Topical Skin Adhesive film to help loosen the bond. Peel off the film; do not pull the skin edges apart.

Disposal

Discard unintentionally opened devices, partially used devices, partially used packages, used devices and packaging according to your facility's policies and procedures concerning biohazardous materials and waste.

Performance / Actions

DERMABOND ADVANCED™ Topical Skin Adhesive is applied topically, with an average setting time of 150 seconds or less, to provide acute incision closure during the normal period of healing, providing a 99% effective microbial barrier for the first 72 hours, then sloughing off within 5-10 days.

Sterility

DERMABOND ADVANCED™ Topical Skin Adhesive is sterilized by dry heat and ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged.

STERILE SINGLE USE ONLY

Storage

Store at or below 30°C (86°F). If the product is exposed outside specified storage range, contact manufacturer. Do not use after expiry date.

How Supplied

DERMABOND ADVANCED™ Topical Skin Adhesive is supplied sterile, in a pre-filled, single-use applicator. The pen style

applicator is composed of a crushable glass ampoule contained within a plastic applicator. The applicator contains the liquid adhesive. The applicator is supplied in a blister package to maintain the device sterility until opened or damaged.

DERMABOND ADVANCED™ Topical Skin Adhesive may be available in boxes of 6 or 12 applicators.

Traceability

The following specific information can be found on the device packaging label: Catalogue number, Batch Code, expiry and manufacturing date, manufacturer name, address and the website and a Unique Device Identification bar code with the Global Trade Item Number information.

Symbols Used on Labeling


	Catalogue number
	Medical Device
	Do not use if package is damaged and consult instructions for use
	Caution
	Do not re-use
	Do not re-sterilize
	Single sterile barrier system with protective packaging inside
	Sterilized using ethylene oxide
	Sterilized using steam or dry heat
	Batch code
	Date of manufacture
	Use-by date
	Unique Device Identifier
	Manufacturer
	Authorized representative in the European Community/European Union

 Rx only

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Packaging unit

www.e-ifu.com
EU: +800 8888 2020
EU: +32 2 4037222

Consult instructions for use or consult electronic instructions for use.

EU: Call paper on demand help desk to get paper copies free of charge within 7 days.

MR Safe

Upper limit of temperature

Contents: Vial

Crush Vial After Opening the Package

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AB: Basılı kopyaları 7 gün içinde ücretsiz olarak talep üzerine almak için danışma hattını arayın.
Ethicon, LLC
Johnson & Johnson Medical GmbH