visocor[®] HM50

Gebrauchsanweisung (DE-2 - DE-37)

Instructions for Use (EN-38 - EN-71)

www.visocor.de

UEBE
Germany est.1890

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Table of contents

A	Intended use	40
B	Safety instructions	
	1. Important instructions for patients	41
	2. Important technical details	43
C	Operating the unit	
	1. Unit description	46
	2. Displays	47
	3. Important instructions for use	47
	4. Initial operation of the device	49
	5. Inserting/replacing batteries, battery icon	49
	6. Activating and setting the date/time	50
	7. Attaching the cuff	51
	8. Body posture when taking measurements	52
	9. Measuring blood pressure	54
	10. Irregular pulses and cardiac arrhythmia	55
	11. Using the memory	56
D	What you should know about blood pressure	
	1. Systolic and diastolic blood pressure values	58
	2. Reasons for measuring different values	58
	3. Why you should measure blood pressure regularly	59
	4. What are normal values for blood pressure?	59

Table of contents

E	Technical information	
	1. Error messages	61
	2. Customer service	63
	3. Technical data	64
	4. Original spare parts and accessories	66
	5. Technical inspection / calibration check	66
	6. Explanation of symbols	67
	7. Disposal	68
F	Maintaining the unit	69
G	Warranty	70

These instructions for use are intended to help the user operate the device safely and efficiently. The unit must be used in accordance with the procedures contained in these instructions for use and must not be used for any other purposes.

It is important that you read all the instructions carefully before using the device. Please pay particular attention to the chapter "Important usage information" on page 47.

A Intended use

It's determined for measuring blood pressure and pulse at the wrist, those of people 15 years of age and over. This unit is recommended for patients with unstable blood pressure for domestic self administration blood pressure control and treatment support.

This unit uses the oscillometric method for measuring blood pressure and pulse rate.

Safety instructions **B**

1. Important instructions for patients

- The device may only be used at the wrist and not on any other extremities.
- Blood pressure measurements on children require specialist knowledge! Please consult your doctor if you wish to measure the blood pressure of a child. Do not under any circumstances use the unit on a baby/infant.
- During pregnancy and in cases of pre-eclampsia, the device may only be used after consulting your doctor.
- Do not under any circumstances place the cuff over a critical area, e.g. a wound, aneurysm, etc. or on an arm with an arteriovenous shunt. Risk of injury! Any supply via an intravascular access point (infusion) or other medical monitoring devices could possibly be interrupted.
- Do not use the unit without first consulting your doctor if you are undergoing dialysis treatment or taking anticoagulants, platelet aggregation inhibitors or steroids. Internal bleeding may be caused under these conditions.

B Safety instructions

- Do not take a blood pressure measurement on the same side as a mastectomy with simultaneous removal of the axilla lymph nodes.
- The display of the pulse frequency is not suitable for checking the frequency of cardiac pacemakers. Cardiac pacemakers and blood pressure monitors do not influence each other regarding their mode of operation.
- The device is not suitable for use in close proximity to electro-surgical instruments or magnetic resonance imaging equipment.
- Do not use the device in explosive environments such as in the vicinity of flammable anaesthetics or in an oxygen chamber.
- Never simultaneously touch the device's live electrical parts and the patient.
- Please wait a few minutes between successive measurements, otherwise the blood flow in the arm is interrupted for too long and injuries could occur.

Safety instructions B

- The device contains small parts that could be swallowed by children. Do not leave the device unattended with children.
- Store the device out of the reach of children.
- Discuss the measured values with your doctor. Do not assess the measurement results yourself. **Never change the dose of medicines prescribed by your doctor.**
- Before carrying out your self-measurement, please pay attention to the chapter "Important usage information" on page 47.

2. Important technical details

- The unit contains sensitive parts and must be protected from extreme temperature fluctuations, humidity, shocks, dust and direct sunlight.
- Do not knock or drop the device.
- The device is not waterproof. Never immerse it in water or any other liquids. If liquid enters the device, this may result in serious damage and may cause malfunctions.

B Safety instructions

- Do not apply pressure to the display. Do not place the device with the display facing downwards.
- The unit must only be operated with the cuff and accessories provided (see “Original spare parts and accessories” on page 66). Using different cuffs or accessories can lead to incorrect measurement results. The warranty will be invalidated if the unit is damaged by third-party accessories!
- A consistently good power supply to the device is necessary for fault-free blood pressure measurements.
 - Only use long-life alkaline batteries (LR03).
 - Always replace all batteries simultaneously during battery replacement. Do not mix new and old batteries or batteries of different types.
 - You will need 2 x 1.5 V AAA/LR03 batteries. Rechargeable batteries have a voltage of only 1.2 V and are thus not suitable.
- If the device is not to be operated for an extended period of time, please remove the batteries. In principle, all batteries can leak.

Safety instructions B

- Using this device in the vicinity of mobile phones, microwaves or other devices with strong electromagnetic fields can cause malfunctions and inaccurate measurements. During use, maintain a distance of at least 30 cm from such devices.
- Never open or modify the device or the cuff (except for replacing the batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- To avoid inaccurate readings, please adhere to the intended operating and storage conditions. See “Technical data” on page 64.
- The inflating and measuring procedure can be interrupted by pressing the Start/Stop button or by removing the cuff. The device then stops the inflation procedure and deflates the cuff.
- If the device malfunctions or is faulty, please refer to the troubleshooting guide starting on page 61, or contact customer service (see page 63).

C Operating the unit

1. Unit description


- 1** Cuff
- 2** Displays
- 3** Start/Stop button
- 4** Memory recall button
- 5** Battery compartment (Removable battery compartment cover)

Operating the unit C

2. Displays

- 1** Storage location or date/time
- 2** SYS = Systole (upper blood pressure value)
- 3** DIA = Diastole (lower blood pressure value)
- 4** PUL 1/min = Pulse calculated pulse frequency per minute
- 5** Pulse signal display or irregular pulse waves
- 6** Battery check display
- 7** Deflating cuff
- 8** Inflating cuff
- 9** Memory identifier

3. Important instructions for use

The measurement results of automated blood pressure monitors can be influenced by the measuring location, posture, previous exertion and general physical condition. Please pay attention to the usage information in order to obtain accurate readings.

C Operating the unit

- Do not consume alcohol, caffeine or smoke for at least one hour before measuring.
- Rest for at least 5 minutes before taking the measurement. Depending on the previous degree of exertion, you may even need to rest for up to an hour.
- Expose the wrist. Under no circumstances may clothes hinder the blood flow to or from the arm, as this affects the blood pressure at the measuring point and can lead to incorrect readings.
- Do not move or talk during measuring.
- Breathe slowly and deeply. Do not hold your breath.
- Check the display for irregular pulse signals. If required, repeat the measurement when conditions have improved.
- In order to exclude side differences and to obtain comparable measurement results, it is important to always take the measurement on the same arm. Ask your doctor which side it is advisable for you to measure.

Operating the unit C

- Blood pressure is not a fixed value. It can go up or down by more than 20 mmHg in patients within a few minutes.

4. Initial operation of the device

Insert the batteries provided into the device.

The device has a date/time function that can be activated as required (see page 50). This function is deactivated on delivery.


When the date/time function is activated, the measurement results are saved along with their date and time. This can make it easier to view the readings at a later point.


5. Inserting/replacing batteries, battery icon

- Remove the cover of the battery compartment on the side of the unit.
- Remove the old batteries from the device and insert the new ones. Please ensure correct polarity (marking in battery compartment).

C Operating the unit

- Close the battery compartment by clipping the battery cover back into the device.

When the “empty battery ” symbol first appears in the display, you are still able to take around 30 measurements. Please replace the batteries within this time frame.

6. Activating and setting the date/time

Activating and deactivating the date/time function

To activate the date/time function, press the Start/Stop button and then immediately, while the full display is visible, the Memory button. This function can be deactivated at any time in the same way.

Setting the date/time

The four-digit number representing the year appears in the display first (Figure 1). Press the Memory button to set the year. Press the Start/Stop but-


Figure 1

Operating the unit C

ton to confirm the selected year and proceed to setting the date (Figure 2).


Figure 2

Set the month with the Memory button and confirm this with the Start/Stop button. Set the day, hours and minutes in the same way (Figure 2, Figure 3).


Figure 3

Changing the date/time

To change the date and time, take the batteries out and wait until the display goes off. Then re-insert the batteries and reset the date and time.


After replacing the batteries, you will have to reset the date and time.

7. Attaching the cuff

The special double-preformed comfort cuff is easy to attach and ensures a perfect fit during measurement. The materials used are tested for skin compatibility according to ISO 10993.

C Operating the unit

- Clear the wrist from all clothing or watches/jewellery since this may disturb the blood flow in the arm as it may impair the blood pressure at the place of measurement.
- Slide the device onto the centre of the inside of the forearm, directly on the wrist. The palm of your hand should be facing upwards. The distance between the wrist and the upper edge of the cuff should be 5–10 mm.
- Fasten the cuff using the hook and loop fastening.

The shape of the cuff is optimized to sit very tightly in order to avoid possible measurement errors.


8. Body posture when taking measurements

The measurements should be performed while seated. In special cases, measurements in a lying position are also possible.

Operating the unit C


Measurement while sitting

Sit at a dining table, lean your back against the backrest of the chair and bend your arm so that your arm is exactly at heart level. Place your feet on the ground and do not cross your legs.

Measurement while lying down

Lay on your back, place the device on your wrist and at the height of your heart, e.g. by supporting it with a pillow. To relax before the measurement, take 2-3 deep breaths.

Please pay attention to the following:

- Do not bend your wrist.
- Keep your hand relaxed, not tightened, do not make a fist.
- It is important to keep your arm calm and not to speak during the measurement.

C Operating the unit

If the cuff is in a higher or lower position than the heart, the measured value might be higher or lower, respectively.

The results for measurements taken in a sitting and lying position may differ slightly.

9. Measuring blood pressure

- Switch on the unit by pressing the Start/Stop button. The full display appears (Figure 1).
- The cuff is inflated to the required pressure. The **▲** symbol flashes in the display (Figure 2).
- The actual measurement process starts with three beeps as the cuff is deflated. The **♥** symbol flashes in the display (Figure 3).
- A long beep denotes the end of measurement. The cuff is automatically deflated.


Figure 1

Figure 2


Figure 3

Operating the unit C

- The determined values of systole, diastole and pulse are shown in the display (Figure 4).
- The device switches itself off after approximately 3 minutes. You can also switch the device off before this by pressing the Start/Stop button.


Figure 4

10. Irregular pulses and cardiac arrhythmia

If the device detects irregular pulses during the measurement, the **♥** symbol is displayed once the measurement is complete.

This can be set off by an irregular heartbeat (arrhythmia), disturbances caused by movement, talking or even breathing in deeply.

If the symbol appears more frequently, this may be caused by an irregular heartbeat and must be discussed with your doctor! Measurement results accompanied by irregular pulses should therefore be viewed as critical and repeated under more favourable conditions.

C Operating the unit

The symbol is stored with the respective measurement.

Caution: Do not confuse the irregular pulse display ▼ with the flashing pulse display ♥ that indicates the heart rate during every measurement (see page 54).

11. Using the memory

Measured results are automatically stored in the memory. The memory can store up to 60 results and the average value.

When more than 60 measured values have been stored, the oldest value (No. 60) is deleted to allow the latest value (No. 1) to be recorded.

Data retrieval

Press the Memory button to call up data. The average value of the stored results is displayed with „A“ (Figure 1).


Figure 1

When the Memory button is pressed again, the value measured last appears (Figure 2). When the

Operating the unit C

date/time function is activated, the measured values are displayed alternately with memory number, date and time.

Press the Memory button repeatedly to call up further measured values.


Figure 2

The stored data remain visible in the display for around 30 seconds. The unit then switches itself off.

Deleting data

To delete individual measured values, press the Memory button repeatedly to call up the desired measured value.

Then press and hold the Memory button again for 8–10 seconds until the measured value flashes several times and then goes out.

To delete the entire memory, retrieve the average value then press and hold the Memory button again for 8–10 seconds until the measured value flashes several times and then goes out.

If no value is stored, there is no display.

D What you should know about blood pressure

1. Systolic and diastolic blood pressure values

The cardiovascular system has the important function of supplying all organs and tissues in the body with sufficient amounts of blood and of transporting metabolites. To do so, the heart contracts and expands at a regular rate of about 60 to 80 times per minute. The pressure of the flowing blood on the artery walls caused by the heart contracting is termed systolic. The pressure in the ensuing relaxation phase, when the heart refills with blood, is termed diastolic. During daily measurement you determine both values.

2. Reasons for measuring different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. Affected by mental, physical and environmental influences, it varies all the time and never remains constant. Reasons for fluctuating blood pressure values can be: moving, speaking, eating, consuming alcohol or nicotine, nervousness, inner tension, emotions, room temperature, recent urination or bowel movements, environmental influences such as movements and sounds, etc. Even changes in the weather and climate can have an effect on your blood pressure. This explains

D What you should know about blood pressure

why values measured at the doctor are often higher than those you obtain at home in your usual environment.

3. Why you should measure blood pressure regularly

Even the time of day has an effect on your blood pressure. The values are at their highest during the day. In the course of the afternoon and in the evening, they drop slightly. They are low while you are sleeping, but rise again relatively quickly once you get up. One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is possible only when measurements are taken regularly. Discuss the measurement values with your doctor.

4. What are normal values for blood pressure?

The World Health Organisation (WHO) established the following limits in mmHg to classify the blood pressure values.

D What you should know about blood pressure

Classification WHO 1999	Systolic pressure = top value mmHg	Diastolic pressure = bottom value mmHg
Stage 3 hypertension	≥ 180	≥ 110
Stage 2 hypertension	160-179	100-109
Stage 1 hypertension	140-159	90-99
High-normal	130-139	85-89
Normal	120-129	80-84
Optimal	< 120	< 80

Depending on age, weight and general health, the blood pressure values can be different.

Only a doctor can determine the right blood pressure range for you and assess whether your blood pressure has reached a dangerous level for you. Discuss your blood pressure values with your doctor. **Never change the dose of medicines prescribed by your doctor!**

Technical information E

1. Error messages

Failure encountered	Possible cause	Corrective action
Display Err - 300	Excessive cuff pressure. Cuff inflated to maximum due to movement of the arm or body.	-Repeat the measurement -Do not move your arm -Do not speak
Display Err - 1	Measurement error; measurement could not be carried out due to movement or speaking.	-Repeat the measurement -Do not move your arm -Do not speak
Display Err - 2	Inflation error, pressure built up too quickly or too slowly.	-Repeat the measurement -Do not move your arm -Do not speak -Check wrist circumference. Call the Service hotline if this occurs repeatedly.
	Cuff is not correctly fitted.	
Display Err - 3	Deflation rate too high or too low.	-Repeat the measurement -Do not move your arm -Do not speak Call the Service hotline if this occurs repeatedly.
	Movement or talking during measurement.	
Display 0 0	Start/Stop button accidentally pressed while battery was being changed.	Switch unit off and then on again with Start/Stop button, if necessary.

E Technical information

Failure encountered	Possible cause	Corrective action
	Irregular heartbeat, disruptive movements, shaking, wobbling (objects), breathing in deeply, etc.	Repeat measurement after 3-5 minutes rest. Possible influences on the measurement results (depending on the severity of the arrhythmia) must be discussed with your doctor!
Measured values/pulse values are extremely high or low.	Too little rest before the measurement. Movement or talking during measurement.	Repeat the measurement after a break of 3-5 minutes. Pay attention to the usage information on page 47.
	Incorrect measurement posture	Repeat the measurement and maintain the correct posture. (Page 52)
Consecutive measurements produce different measured values (readings).	Blood pressure is not a fixed quantity. It may vary in the upward or downward direction within a few minutes by more than 20 mmHg.	In order to obtain comparable results, always carry out the measurements under the same conditions.
Measurement was interrupted, the cuff deflated to 0 and the measurement restarted.	Device detects a fault or that the diastolic value is too low.	Device shows no fault. Repeat measurement
	Movement during measurement	If necessary, interrupt measurement, measure again after a break of 5 minutes.

Technical information E

Failure encountered	Possible cause	Corrective action
No display after unit is switched on	Batteries inserted incorrectly.	Check position of batteries.
	Batteries flat.	Change batteries.
	Battery contacts dirty.	Clean battery contacts with a dry cloth.
Nothing appears in the display after pressing the Memory button	No measured values have been stored.	Repeat measurement
Measurement interrupted	Batteries flat.	Change batteries.
Date/time not shown in display	Function accidentally deactivated	Activate date and time as described on page 50.

2. Customer service

Device repairs may only be performed by the manufacturer or an expressly authorised body. Please contact:

E Technical information

UEBE Medical GmbH
 Zum Ottersberg 9
 97877 Wertheim, Germany
 Tel.: +49 (0) 9342/924040
 Fax.: +49 (0) 9342/924080
 Email: info@uebe.com
 Website: www.uebe.com

3. Technical data

Model:	visocor® HM50
Size:	approx. W = 70 mm x H = 70 mm x D = 80 mm with cuff
Weight:	121 g with cuff, not including batteries
Display:	LCD display (liquid crystal display) for measured values and check displays
Memory:	60 measured values (stored automatically) and average value (A)
Measurement procedure:	Oscillometric determination of systole, diastole and pulse
Reference procedure of clinical testing:	Auscultatory measurement
Inflating pressure:	Automatically controlled
Pressure display range:	0-300 mmHg

Technical information E

Measurement range:	Systolic: 50-250 mmHg Diastolic: 40-180 mmHg Pulse measurement: 40-160 pulses/min
Measurement precision:	Blood pressure measurement (cuff pressure): ± 3 mmHg Pulse rate: ± 5%
Power supply:	Battery type: 2 x 1.5 V, AAA, LR 03, durability: over 700 measurements in 2 years
Cuff:	Cuff type H5, wrist circumference 12.5 cm - 21.5 cm
Operating conditions:	Ambient temperature 10°C to 40°C, relative humidity 15% to 85% non-condensing, air pressure 700 to 1060 hPa
Storage and transportation conditions:	Ambient temperature -5 °C to 50 °C Relative humidity up to 85%, non-condensing
Inflation and deflation:	Regulated automatically
Protection class IP:	IP 20: Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.
Protection against electric shock:	Internal power supply, applied part type BF (cuff)
Expected service life:	5 years
Operating mode:	Continuous operation
Classification:	Internally powered ME Equipment

E Technical information

4. Original spare parts and accessories

The following original spare parts and accessories are available from specialist dealers:

- Cuff type H5 (for wrist circumference 12.5-21.5 cm)
Part no. 2205001, PZN-09101582

Subject to technical modifications.

5. Technical inspection / calibration check

Generally, it is recommended that a metrological inspection is performed every two years. However, professional users in Germany are obligated to do so in accordance with the "Medical Products Operator Ordinance" (Medizinprodukte-Betreiberverordnung). This can be performed either by UEBE Medical GmbH, an authority responsible for metrology or an authorised maintenance service. Please refer to your national regulations.

Upon request, responsible authorities or authorised maintenance services receive a „Test instruction for metrological inspection“ from the manufacturer.

Technical information E

Please only submit the device for metrological inspection together with the cuff and the instructions for use.

Important: No modifications, e.g. opening the device (except to replace the batteries), may be made to this device without the manufacturer's permission.

6. Explanation of symbols

This product complies with Council Directive 93/42/EEC of 5 September 2007 concerning medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Degree of protection against electric shock TYPE BF

Consult the instructions for use

Storage and transportation conditions

Ambient temperature -5 to +50 °C

Relative humidity up to 85 %, non-condensing

E Technical information

Keep dry

SN

Serial number

REF

Reference number = item number

IP20

Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.

Manufacturer

7. Disposal

Batteries and technical appliances must not be disposed of with domestic waste, but should be handed in at the appropriate collection and disposal points.

Maintaining the unit F

- Keep the device clean. Check for cleanliness after use. Please use a soft, dry cloth for cleaning. Do not use benzene, thinners or other strong solvents.
- The cuff may absorb sweat and other fluids. Inspect the cuff for marks or discolourations after each use. Please use a synthetic detergent for cleaning and carefully wipe the surface. Do not scrub or machine-wash. Air-dry thoroughly.
- Make sure when storing the unit that no heavy objects are left lying on the unit.
- Do not apply force when folding or bending the cuff.
- If the unit has been stored at temperatures below 10 °C or above 40 °C, let it stand at room temperature for at least 2 hours before use.
- If the device is not to be operated for an extended period of time, please remove the batteries. Store the batteries out of the reach of children.

G Warranty

The device has been manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we provide warranty in accordance with the following terms and conditions:

During the warranty period of 2 years from the date of purchase we reserve the right either to repair any such defect at our expense or to supply a perfect replacement unit. The cost of returning the unit to our factory shall be borne by the sender.

Excluded from the warranty are parts subject to normal wear and tear as well as damage caused by non-compliance with the instructions for use, improper handling (e.g. unsuitable power sources, breakages, leaking batteries) and/or disassembly of the unit by the purchaser. Furthermore, no claims for damages against us are substantiated by the warranty.

Warranty claims can only be advanced in the warranty period and by presenting proof of purchase. In the event of a warranty claim, the unit must be sent to the following address together with the proof of purchase and a description of the complaint:

Warranty **G**

UEBE Medical GmbH
Service-Center
Zum Ottersberg 9
97877 Wertheim, Germany

The cost of returning the device to our factory shall be borne by the sender. Complaints that are sent back without prepayment will not be accepted by UEBE.

The statutory claims and rights of the buyer against the seller (for example, defect claims, producer liability) are not restricted by this warranty.

Please note: In the case of a warranty claim, please make sure to enclose the proof of purchase.

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© Copyright 2018 UEBE Medical GmbH