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### DISPOSITIF MEDICAL

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Solution ophtalmologique 15 ml spray

## EN eye.care//comfort spray Eye spray with liposomes and hyaluronic acid 0,1%.

Preservative free lubricating, moisturizing and refreshing eye and eyelid spray with Liposomes and Hyaluronic acid 0,10% for replenishing and stabilizing the tear film. Helps to improve the symptoms of dry eyes. Provides relief to the eye, eyelid and periorcular area in case of sensitive, dry and dehydrated eyes. Particularly indicated for contact lens wearers. The liposome eye spray is based on a phospholipid complex which encapsulates Hyaluronic acid. Liposomes are a typical example of colloidal vector, mainly made up of phospholipids, the same molecules which make up the cellular membrane matrix, thus determining a high biocompatibility. Its spray technology ensures a better diffusion and hydration than traditional drops.

The product remains sterile after first use.

### INDICATION

The liposome eye spray carries out a lubricating and moisturizing action, useful in cases of dry eye syndrome and of alterations of the lipidic layer of the tear film, particularly for contact lens wearers, thanks to transporting the Hyaluronic acid and phospholipids to the eye. Provides comfort to eyelids and relieves the sensation of dry eye, foreign body, burning and itching eye and eyelids. Hyaluronic acid, in the form of Sodium hyaluronate, is a biological polysaccharide with the same molecular structure of the intraocular fluids and, thanks to its viscoelasticity and hygroscopic characteristics, it can significantly improve the tear film stability. The vehiculation of Hyaluronic acid at eye level enhances hydration and reduces the risk of dry eye sensation. As a mucomimetic, Hyaluronic acid has the ability to simulate some of the properties of the natural tear mucin to correct the instability of the tear film. Liposomes (phospholipid complex) reconstitute and reintegrate the lipid layer of the tear film, that is the protective barrier which controls its evaporation. The combined action of phospholipids and Hyaluronic acid is therefore stabilizing and has the purpose of maintaining and restoring the stability of the tear film, ensuring moistening of the exposed cornea-conjunctival surface and a uniform positioning of the covering tear film.

### INSTRUCTIONS FOR USE

First use and learning phase: the product is administered with a no gas spray system; therefore, it may be necessary to press the button repeatedly to obtain the first delivery. On this occasion, it is advisable to direct the spray towards the palm of your hand in order to check the proper degree of pressure to be applied to obtain optimal atomization. Pressing the button indecisively will not have an ideal effect. Thoroughly wash your hands. Shake before use. Hold the bottle tightly so that your thumb (or forefinger) is firmly on the button, being careful not to cover the outlet hole of the dispenser. Exert decisive pressure from a distance of about 10-15 cm from the eye in a slightly lateral position on the closed eye, so that the spray is evenly distributed not only on the eye, but also on the eyelid area, on the area around the eyes and on the eyebrows. It is recommended to spray twice at each application, once or more times a day. Do not use with make-up which isn't waterproof, the liposomes could vehiculate component of makeup into the eyes.

### WARNINGS

In case of eye irritation stop using and contact a specialist. Do not use in case of allergy or hypersensitivity to any of the product components. Wash your hands before each application. Keep away from heat sources, store in a cool, dry place at a temperature between 5 and 35°C. Do not use the product if the bottle appears to be damaged. Close the bottle firmly after each use. Keep out of the reach of children. Do not swallow. Do not use the product after the expiry date. Use the product within 90 days of first use. Dispose of the packaging according to the local laws.

### COMPOSITION

Phospholipidic complex, Hyaluronic acid (as Sodium hyaluronate) 0.10%. Boric acid, Sodium tetraborate. Isotonic solution buffered to pH 7.2. Contains soy lecithin. Preservative free.

### MEDICAL DEVICE

### FORM AND CONTENT

Ophthalmic solution in a 15 ml spray bottle.

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