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Text Rev.06 21.06.22

www.schalcon.com

UK SINGLE-PHASE PEROXIDE SYSTEM WITH PLATINUM CATALYST FOR CONTACT LENS CARE PROPERTIES

Single-phase peroxide is a system suitable for regular contact lens care. Regular use of the product avoids the use of complementary systems for contact lens maintenance, such as detergents and/or enzymatic tablets. As the product is a hydrogen peroxide system, it disinfects and cleans contact lenses removing proteins, fats and all organic residues that adhere to contact lenses when worn. Furthermore, once the hydrogen peroxide is eliminated by catalysis with platinum, the product hydrates the lenses, making them ready for use. Daily use of the Single-phase Peroxide System with platinum catalyst allows using contact lenses which are always disinfected, perfectly clean and hydrated. The product is a peroxide system activated by the platinum catalyst that combines its exceptional antiseptic and ultra-cleaning abilities with the destructive action of nascent oxygen on bacteria, fungi, viruses and organic residues which may adhere to contact lenses when worn. The gaseous bubbles, generated when the solution comes into contact with the platinum catalyst, are a sign that the product has started its disinfectant and ultra-cleaning action.

How it works

When used in combination with the platinum catalyst disc (pre-inserted in the lens case), the peroxide solution develops a disinfectant action on the contact lenses. After 6 hours the neutralisation process is complete and the peroxide solution is completely transformed into a preservative-free saline solution. The solution must always be used in combination with the lens case and catalyst disc contained in this package. Failure to use the catalyst disc may cause burning and/or sharp eye pain.

Taking care of your eyes

A proper lens disinfection is a fundamental step in the lens care regimen. The single-phase hydrogen peroxide solution is a simple method for disinfecting and storing contact lenses, particularly suitable for regular users. It is very important to remember that contact lenses need proper maintenance to help keep your eyes healthy. For this reason, the instructions for use must be followed scrupulously.

CONTENT

One bottle of sterile, buffered solution, containing 3% hydrogen peroxide.

A lens case with pre-inserted catalyst disc.

A leaflet.

INDICATIONS

The single-phase peroxide system with platinum catalyst is suitable for care of all contact lens types: it disinfects, cleans, hydrates and removes proteins, fats

and all organic residues that adhere to contact lenses.

INSTRUCTIONS FOR USE

- Always wash and dry your hands before handling the lenses.
- Unscrew the cap of the special lens case supplied and check that the catalyst disc is present.
- Put the left contact lens in the basket marked with the symbol with the letter L. Close the basket taking care not to damage the lens.
- Repeat the same operation with the right lens, placing it in the other basket marked with the letter R.
- Fill the lens case with peroxide solution up to the reference line. Do not fill it above or below the indicated reference line. First place the lenses in the baskets, so that they are ready to be inserted into the lens container with the solution. This will prevent the neutralisation process from starting without the lenses.
- Place the support with the baskets containing the lenses into the case and screw it tightly into place.
- The catalyst disc will begin to produce bubbles, thus showing that the neutralisation process has begun.
- Leave the lenses in the peroxide solution for at least 6 hours.
- At the end of the treatment, gently shake the lens case and then unscrew the lid and remove the lenses first from one basket and then from the other.
- Always rinse the lenses with saline solution before wearing them.
- Discard any solution left in the lens case immediately. Leave the lens case open to air dry.

Your contact lens specialist may recommend cleaning your lenses daily with a specific cleaner. In this case it is necessary to rinse and completely remove all traces of the detergent from the lenses with the saline solution, before immersing them in the peroxide solution.

For some patients, it may be necessary to remove protein deposits by using enzymatic tablets. In this specific case, consult the contact lens specialist and follow his/her instructions.

PRECAUTIONS FOR USE

- Always use a new lens case with the pre-installed catalyst disk contained in the package. Discard the bottle and the lens holder 3 months after the first use or in any case when the solution is finished.
- Follow the maintenance procedures recommended by the contact lens specialist. Failure to follow these instructions can cause serious damage to eyes.
- As soon as the solution comes into contact with the platinum catalyst, it will release some gas bubbles. This reaction is the confirmation that the catalyst is activating the Hydrogen peroxide. Should this not occur, lightly shake the container and wait for the reaction to start.
- If the lenses have been immersed in peroxide solution for more than 24 hours, disinfect them again, replacing the solution in the lens case with fresh peroxide solution, following the instructions in this package leaflet.
- Do not fill the lens case with peroxide solution above or below the reference line indicated.
- Fill the lens case with peroxide solution only after the lenses have been placed in the baskets. This will prevent the neutralisation process from starting without the lenses.

COMPOSITION

100ml of solution contains:

Hydrogen peroxide 3.00 g
Potassium chloride 0.06 g
Pluronic RPE 1740 0.008 g
Stabiliser 0.015 g

Quantum satis buffered isotonic solution at physiologic pH.

Catalyst included in the lens holder: ABS (Acrylonitrile Butadiene Styrene) coated with Platinum.

WARNINGS ⚠

- Do not use peroxide solution directly in contact with eyes.
- Always rinse the contact lenses with saline solution before wearing them.
- Always wait at least 6 hours before wearing the contact lenses.
- Use the peroxide solution only with the lens case provided, making sure that the platinum catalyst is always inside.
- Avoid use in temperatures below 15 °C.
- If the non-neutralised peroxide solution comes into contact with the eyes, rinse immediately with running water.
- Only use fresh peroxide solution for each disinfection cycle and discard the remaining solution left in the lens case after each use.
- Do not use the product after the expiry date shown on the package.
- Use the product within 3 months of first opening.
- Do not swallow. Do not freeze, keep away from heat sources. KEEP OUT OF REACH OF CHILDREN. Tightly close the bottle after use.
- Do not use in case of known hypersensitivity to one of the components.

IMPORTANT

If you experience irritation or persistent or increasing excessive tearing, or if your vision is impaired, stop wearing your contact lenses and consult your ophthalmologist.
PRESERVATIVE FREE
Medical Device

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