## LadySystem Vaginalkonen

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### Produzent:

LADYSYSTEMS INTERNATIONAL LIMITED  
Waldfield House  
32 High Street  
Pinner, HAA5 5PW  
UK

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## LadySystem Vaginal Konen

### The importance of the musculature of the pelvic floor

The musculature of the pelvic floor acts as a cradle for the bladder, the uterus and the rectum/anterior intestine and keeps these organs in their correct position (picture 1). In order for these organs to function correctly, the musculature of the pelvic floor must be healthy and strong.

Picture 1: Pelvic floor in normal condition

### The factors which cause weakness of the pelvic floor

There are different reasons for weakness of the pelvic floor. Pregnancy and childbirth are the main causes. During pregnancy the increased intra-abdominal weight places a large additional strain on the pelvic floor stretching it, and this reaches a maximum during birth. The hormonal changes of the menopause can result in a loss of flexibility in the pelvic floor. Other physical stress such as sports which involve jumping or running (e.g. aerobic, tennis) or playing wind instruments can also stretch and weaken the musculature of the pelvic floor.

### Effects which pelvic floor weakness can provoke

Stress incontinence is a very common consequence of weakness of the pelvic floor. It is the leakage of a variable quantity of urine after a sudden increase in intra-abdominal pressure, which takes place for

example on coughing, sneezing, laughing, walking or simple jumping.

### Also descent of the uterus (uterine prolapse) or bladder could be caused by a weakness of pelvic floor.

### VAGINAL CONES

Vaginal cones have been developed as a simple and effective method for the non-surgical treatment of weak pelvic floor muscles and for the prevention of the conditions associated with this problem. LadySystem vaginal cone therapy consists of five anatomically shaped cones, each of slightly different weight, which are used by inserting the appropriate one in the vagina for 15 minutes twice a day. They work while the user is standing or walking, allowing the continuation of normal daily tasks. An improvement in the pelvic floor tone is noticeable within 2 or 3 weeks and the majority of women recover full tone after between 6 weeks and 3 months.

### Indications for vaginal cones

- Treatment of Stress Incontinence**  
If already present to greater or lesser degree, this is indicative of a weak pelvic floor and it is advisable to seek treatment before the condition deteriorates further.
- Auxiliary therapy in surgical procedures**  
Preoperative rehabilitation of the pelvic floor muscles in patients awaiting surgery, in many cases this can obviate the need for surgery. Used as a post-operative therapy, vaginal cones can improve the result of surgery and help to prevent relapse.

### Re-education of the musculature post-partum

The effectiveness of cone therapy in restoring pelvic floor tone has led to this technique being routinely recommended as a preventive therapy for all mothers.

### Evaluation of the condition of the pelvic floor and maintaining pelvic floor tone

Cones can serve as an objective measure


of pelvic floor strength. It is important to carry out periodic evaluation of the pelvic floor, and it is recommended that this is done every six months. The exercises can be re-started at any time if the pelvic floor has weakened again.

### INSTRUCTIONS FOR USE

The set of LadySystem consists of five tampon-shaped cones with different weights (4.9 g; 20.5 g; 32 g; 44 g; 55 g). The cones are for personal use only and should not be used by any other person, in order to avoid passing on infections. Before and after each use, use soap and water to wash the cones which is to be used. After use, dry the cone with a clean towel before storing in the holder. It is not necessary to sterilise the cones, but if you wish to do so, use a chemical solution of the type used for sterilising baby bottles. They can be damaged if in contact with hot water or alcohol.

Unless directed by your doctor, it is not advisable to use the cones if you have a vaginal infection; during pregnancy, while an episiotomy wound is healing or is causing any discomfort (if in doubt consult your doctor).

- Begin with the cone No. 1, the lightest one. Insert it smoothly into the vagina just like a tampon. You can moisten it with water to ease its insertion. The cone should be completely inside and nearly vertical in the vagina, with only the nylon cord outside (picture 2).

Picture 2: Correct insertion